12 new publications - Vapourtec Winter 2014 newsletter



A warm welcome to the Vapourtec Winter 2014 e-newsletter.

A Merry Christmas and a very happy 2015 from all of the team at Vapourtec!

Contents

Product News

• Introducing the R2-S new Suspension / light slurry pump module

Latest News

- Continuous flow of publications tops 150 for Vapourtec systems
- Vapourtec sees growth in India
- Flow v batch Photochemistry a like-for-like comparison

Events

• Meetings, conferences or exhibitions where you can see Vapourtec products and talk to Vapourtec staff.

Applications

- Synthesis of Artemisinin via the Photooxidation of Dihydroartemisinic Acid
- Preparation of Silver Nanoparticles under Continuous Flow Conditions

Publications

• 12 new publications in this issue

Product News

Introducing the new R2-S Suspension / light slurry pump module

Vapourtec have announced a new addition to the range of pump modules available for the R-SeriesTM system. In addition to the standard, acid resistant and high pressure there is now a special suspension / slurry capable version.

The new module has all the same control and data logging functionality as the standard pump with some additional pumping capabilities...

MORE DETAILS



Latest News

A continuous flow of published breakthroughs for Vapourtec

Vapourtec have recently reached the impressive milestone of having been cited in 150 peer review publications.

By way of a comparison Vapourtec's two closest competitors in flow chemistry technology have been mentioned in 51 and 40 publications respectively. READ MORE



Vapourtec sees flow chemistry growth in India

Vapourtec is establishing a burgeoning market for its flow chemistry systems across India with the continued development of a distribution partnership Mumbai-based **Pi-Process** with Intensification. Pi-Process Intensification have been acting as distributors for Vapourtec's R-Series systems for the past three years and will now also be distributing their E-Series system. Dr Reddy's Labs is one of a number of key customers. **READ MORE**



Flow v Batch Photochemistry - a like-for-like comparison

There are a number of key advantages of a continuous flow approach to photochemistry over a traditional batch method including consistent light penetration, controlled exposure times, precise temperature control and easy scalability as well as the removal of photochemical products from the irradiated area.



READ MORE

Events

Events where you can see Vapourtec systems in action:

15th – 16th December 2014 RSC Macrocyclic and Supramolecular Chemistry Meeting Norwich, UK Click <u>here</u> for more details

22nd – 23rd January 2015 Flow Chemistry India 2015 Mumbai, India Click <u>here</u> for more details

17th – 18th February 2015 Flow Chemistry Conference Europe Berlin, Germany Click <u>here</u> for more details

22nd – 26th March 2015 ACS Spring Meeting Denver, USA Click <u>here</u> for more details

1st April 2015 Reagentless Synthesis SCI, London, UK Click <u>here</u> for more details

Applications

H

Synthesis of Artemisinin via the Photooxidation of Dihydroartemisinic Acid

	1) O _{2, hv} , PhCH ₃					
	15 °C, 5 mins 2) 25 °C, 15 mins					
Ĥ	Ī	Ĥ	The second secon	_	H	10000



Scheme 1: Synthesis of artemisinin from dihydroartemisinic acid

Application Note 39: This Application Note shows the versatility of the Vapourtec easy-PhotoChem system. Artemisinin is synthesized directly from dihydroartemisinic acid via the continuous photooxidation route first published by Seeberger et al in 2012. In this example the easy-PhotoChem system is used to, pump liquid reagents, meter oxygen and provide the light source for the photooxidation. Artemisinin was produced at an impressive rate of 1.6 grams / hour.

Preparation of Silver Nanoparticles under Continuous Flow Conditions



Application Note 40: This Application Note describes the controlled formation of silver nanoparticles in tubular reactors using the Vapourtec E-Series. Control of particle size is shown over the size range 10 to 60 nm. Two classes of nanoparticle are reported.

Click here to go to the Application Notes page on the Vapourtec website

Publications

Chemical Assembly Systems: Layered Control for Divergent, Continuous, Multistep Syntheses of Active Pharmaceutical Ingredients



While continuous chemical processes have attracted both academic and industrial interest, virtually all active pharmaceutical ingredients (APIs) are still produced by using multiple distinct batch processes. To date, methods for the divergent multistep continuous production of customizable small molecules are not available. A chemical assembly system was developed, in which flow-reaction modules are linked together in an interchangeable fashion to give access to a wide breadth of chemical space. Control at three different levels—choice of starting material, reagent, or order of reaction modules—enables the synthesis of five APIs that represent three different structural classes (γ -amino acids), including the blockbuster drugs Lyrica and Gabapentin, in good overall yields (49–75%).

Click here to go straight to the publication

Continuous Synthesis of Organozinc Halides Coupled to Negishi Reactions Nerea Alonso^{2,3},

L. Zane Miller¹,

Juan de M. Muñoz²,

Jesus Alcázar^{2,*}

D. Tyler McQuade^{1,*}

¹Department of Chemistry and Biochemistry, Florida State University, USA ²Janssen Research and Development, Janssen-Cilag, Toledo, Spain ³Facultad de Química, Universidad de Castilla-La Mancha, Spain

The Negishi cross-coupling is a powerful CC bond forming reaction. The method is less commonly used relative to other cross-coupling methods in part due to lack of availability of organozinc species. While organozinc species can be prepared, problems with reproducibility and handling of these sensitive species can complicate these reactions. Herein, we describe the continuous formation, using an activated packed-bed of metallic zinc, and subsequent use of organozinc halides. We demonstrate that a single column of zinc can provide excellent yields of organozinc halides and that they can be used downstream in subsequent Negishi cross-couplings. The preparation of the zinc column and the scope of the reaction are discussed.

Click here to go straight to the publication

Continuous Flow Magnesiation of Functionalized Heterocycles and

Acrylates with TMPMgCl·LiCl.



A flow procedure for the metalation of functionalized heterocycles (pyridines, pyrimidines,



thiophenes, and thiazoles) and various acrylates using the strong, non-nucleophilic base TMPMgCI-LiCI is reported. The flow conditions allow the magnesiations to be performed under more convenient conditions than the comparable batch reactions, which often require cryogenic temperatures and long reaction times. Moreover, the flow reactions are directly scalable without further optimization. Metalation under flow conditions also allows magnesiations that did not produce the desired products under batch conditions, such as the magnesiation of sensitive acrylic derivatives. The magnesiated species are subsequently quenched with various electrophiles, thereby introducing a broad range of functionalities.

Click here to go straight to the publication

Continuous Synthesis of Artemisinin-Derived Medicines

Kerry Gilmore,^a

Daniel Kopetzki,^a

Ju Weon Lee,^b

Zoltan Horvath,^b

D. Tyler McQuade,^a

Andreas Seidel-Morgenstern,^{b,c}

Peter H. Seeberger^{a,d}

^a Max-Planck-Institute of Colloids and Interfaces, Department of Biomolecular Systems, Germany

^b Max-Planck-Institute for Dynamics of Complex Technical Systems, Germany

^c Otto-von-Guericke-University, Chair for Chemical Process Technology, Germany

^d Freie Universität Berlin, Institute of Chemistry and Biochemistry, Berlin, Germany

Described is a continuous, divergent synthesis system which is coupled to continuous purification and is capable of producing four anti-malarial APIs. The system is comprised of three linked reaction modules for photooxidation/cyclization, reduction, and derivatization. A fourth module couples the crude reaction stream with continuous purification to yield pure API.

Click here to go straight to the publication



Continuous flow macrocyclization at high concentrations: synthesis of macrocyclic lipids

Anne-Catherine Bédard, Sophie Régnier, Shawn K. Collins



Département de Chimie, Centre for Green Chemistry and Catalysis, Université de Montréal, Montréal, Canada

A phase separation/continuous flow macrocyclization protocol eliminates the need for high-dilution conditions and can be used to prepare gram quantities of biologically relevant macrocyclic lipid structures. The method presents several green advantages towards macrocycle synthesis: (1) the prevention of unwanted oligomers and waste, (2) a reduction in the large quantities of toxic, volatile organic solvents and (3) the use of PEG as an environmentally benign reaction media. Macrocycles could be synthesized in high yields (up to 99%) in short reaction times (1.5 h) and on gram scales without the need to alter the reaction conditions.

Click here to go straight to the publication

Versatile, High Quality and Scalable Continuous Flow Production of

Metal-Organic Frameworks



CSIRO Materials Science and Engineering, Australia

Further deployment of Metal-Organic Frameworks in applied settings requires their ready preparation at scale. Expansion of typical batch processes can lead to unsuccessful or low

quality synthesis for some systems. Here we report how continuous flow chemistry can be adapted as a versatile route to a range of MOFs, by emulating conditions of lab-scale batch synthesis. This delivers ready synthesis of three different MOFs, with surface areas that closely match theoretical maxima, with production rates of 60 g/h at extremely high space-time yields.

Click here to go straight to the publication

Still reading someone else's copy of the newsletter?



If you like reading the newsletter but tend to get it forwarded by a colleague, why not sign up for your own copy? It will only take a minute, and your email address won't be used for anything else, ever. Each newsletter contains an "unsubscribe" link in case you should change your mind.

Click here to signup for your own copy

First Example of a Continuous-Flow Carbonylation Reaction Using Aryl Formates as CO Precursors

Nerea Alonso^{1, 3},

Juan de Muñoz¹,

Brecht Egle²,

Johannes L. Vrijdag²,

Wim M. De Borggraeve²,

Antonio de la Hoz³,

Angel Díaz-Ortiz³,

Jesús Alcázar¹

¹Janssen Research and Development, Janssen-Cilag Department of Medicinal Chemistry S.A., Toledo, Spain

²Molecular Design and Synthesis University of Leuven, Department of Chemistry, Heverlee Belgium

³Universidad de Castilla-La Mancha Facultad de Ciencias y Tecnologías Químicas Real,



The first continuous flow carbonylation reaction using aryl formates as CO precursor is reported. The reaction is practical, scalable and high yielding. The use of a flow protocol safely allows expanding the scope to activated chlorides, nitrogen heterocycles and to the selective introduction of an ester group in dihalo-derivatives. Further selective reduction of the ester formed to an aldehyde in flow is also described.

Click here to go straight to the publication

Multistep Flow Synthesis of 5-Amino-2-aryl-2H-[1,2,3]-triazole-4-

carbonitriles

Dr. Jérôme Jacq,

Dr. Patrick Pasau*



UCB Biopharma, Avenue de l'Industrie, 1420 Braine l'Alleud (Belgium)

1,2,3-Triazole has become one of the most important heterocycles in contemporary medicinal chemistry. The development of the copper-catalyzed Huisgen cycloaddition has allowed the efficient synthesis of 1-substituted 1,2,3-triazoles. However, only a few methods are available for the selective preparation of 2-substituted 1,2,3-triazole isomers. In this context, we decided to develop an efficient flow synthesis for the preparation of various 2-aryl-1,2,3-triazoles. Our strategy involves a three-step synthesis under continuous-flow conditions that starts from the diazotization of anilines and subsequent reaction with malononitrile, followed by nucleophilic addition of amines, and finally employs a catalytic copper(II) cyclization. Potential safety hazards associated with the formation of reactive diazonium species have been addressed by inline quenching. The use of flow equipment allows reliable scale up processes with precise control of the reaction conditions. Synthesis of 2-substituted 1,2,3-triazoles has been achieved in good yields with excellent selectivities, thus providing a wide range of 1,2,3-triazoles.

The direct α -C(sp³)–H functionalisation of N-aryl tetrahydroisoquinolines via an iron-catalysed aerobic nitro-Mannich reaction and continuous flow processing.

Martin Brzozowski,		FeCl ₂ (10 mol%) RCH ₂ NO ₂ (5 equiv.)			
Jose A. Forni,	\cap		\bigcirc		
G. Paul Savage	`` ^N `Ar		````Ar		
Anastasios Polyzos*	п	O₂ (7 bar), MeOH 90 °C 0.167 mL min ⁻¹	R NO ₂		
	Teflon AF-2400 Tube-in-Tube				

CSIRO Manufacturing Flagship, Bayview Avenue, Clayton, Australia

An efficient nitro-Mannich type direct α -C(sp³)–H functionalisation of *N*-aryl-1,2,3,4-tetrahydroisoquinolines catalysed by simple iron salts in combination with O₂ as the terminal oxidant is described. The use of a Teflon AF-2400 membrane Tube-in-Tube reactor under continuous flow conditions allowed for considerable process intensification to be achieved relative to previous batch methods.

Click here to go straight to the publication

A Continuous-Flow Approach to 3,3,3-Trifluoromethylpropenes: Bringing Together Grignard Addition, Peterson Elimination, Inline

Extraction, and Solvent Switching



[†] Department of Chemistry, University of Connecticut, United States

[‡] Department of Community Medicine & Health Care, University of Connecticut Health Center, United States

A continuous-flow approach to the synthesis of 3,3,3-trifluoromethylpropenes involving Grignard addition of (trimethylsilyl)methylmagnesium chloride to a trifluoromethyl ketone followed by dehydrative desilylation of the α -trifluoromethyl- β -hydroxysilyl alcohol using trimethylsilyl trifluoromethanesulfonate is reported. An inline aqueous/organic extraction and a concomitant solvent switch were key to the success of the methodology. Transition from batch to continuous flow conditions allows for higher yields, shorter reaction times, and facile scale out.

Click here to go straight to the publication

Continuous-Flow Oxidative Cyanation of Primary and Secondary Amines

O2, TPP, THF

LED 420 nm TMSCN, r.t. $|_{Ph} \sim_{N} \sim_{Ph} | \rightarrow Ph^{\perp}$

99%

Using Singlet Oxygen

Dmitry B. Ushakov,

Kerry Gilmore,

Daniel Kopetzki, D.

Tyler McQuade,

Peter H. Seeberger

¹Department für Biomolekulare Systeme, Max-Planck-Institut für Kolloid- und Grenzflächenforschung, Potsdam, Germany

²Institut für Chemie und Biochemie, Freie Universität Berlin, Berlin, Germany

³Department of Chemistry and Biochemistry, Florida State University, Tallahassee, USA

Primary and secondary amines can be rapidly and quantitatively oxidized to the corresponding imines by singlet oxygen. This reactive form of oxygen was produced using a variable-temperature continuous-flow LED-photoreactor with a catalytic amount of tetraphenylporphyrin as the sensitizer. α -Aminonitriles were obtained in good to excellent yields when trimethylsilyl cyanide served as an in situ imine trap. At 25 °C, primary amines were found to undergo oxidative coupling prior to cyanide addition and yielded secondary α -aminonitriles. Primary α -aminonitriles were synthesized from the corresponding primary amines for the first time, by an oxidative Strecker reaction at -50 °C. This atom-economic and protecting-group-free pathway provides a route to racemic amino acids, which was exemplified by the synthesis of *tert*-leucine hydrochloride from neopentylamine.

Click here to go straight to the publication

Facilitating Biomimetic Syntheses of Borrerine Derived Alkaloids by

Means of Flow-Chemical Methods.



^A Department of Chemistry, University of Cambridge, Lensfield Road, Cambridge CB2 1EW, UK.

^B Corresponding author.

Flow chemistry is widely used nowadays in synthetic chemistry and has increasingly been applied to complex natural product synthesis. However, to date flow chemistry has not found a place in the area of biomimetic synthesis. Here we show the syntheses of borrerine derived alkaloids, indicating that we can use biomimetic principles in flow to prepare complex architectures in a single step.

Click here to go straight to the publication



Thanks for reading.

See you again in 2015!

f Share Image: Share Image: Share Image: Share Image: Share
f 💟 in
Facebook Twitter LinkedIn

Copyright © 2014 Vapourtec Ltd, All rights reserved. unsubscribe from this list update subscription preferences