

Introducing the R2-S – The R-Series slurry capable pump module

In addition to the standard, acid resistant and high pressure pump modules for the R-Series, Vapourtec, has now launched a suspension and light slurry version; **R2-S**.

The **R2-S** features Vapourtec's high pressure precision peristaltic **V-3** pump - capable of pumping at up to 10bar, with flow rates from 0.02 to 10ml/min.

Robust:


Simple, reliable and smooth pumping across a wide range of flow rates

Capable of pumping:

- Strong acids
- Organometallics
- Suspensions & light slurries
- Gases

Easy to use:

- Self priming pumps
- Adjustable acid compatible back pressure regulator
- Can pump two consecutive solvents even if they are not miscible




The Vapourtec **V-3** features advanced control for smooth output.

The tube materials are specially developed in collaboration with a major global fluoropolymer manufacturer, resulting in a pump which can generate high pressures and is compatible with synthesis solvents, acids and bases.

How does the Vapourtec V-3 work ?

A standard peristaltic pump runs at a steady rotational speed whereas the V-3 does not.

The exact rotational speed profile is complex and is a patented function of the flow rate, rotor position and pressure. The result of this advanced development is total simplicity for the user.



There are two configurations of the new slurry capable pump module:

- **R2-S** without sample injection loops and the
- **R2-S+** with sample injection loops

The system can be controlled manually with the built in interface or can be used with FlowCommander software, which allows automation of experiments, logging, charting, reporting, saving and sharing of reaction setup and data.

The **R2-S** can be integrated with other equipment to provide versatile automated flow chemistry including Autosampler, UV and FTIR.




The **R2-S** can be added to an existing R-Series system to extend capabilities or can be used as the primary pump module for a system. This will still allow expansion of your flow chemistry capability as your needs develop.

None of your investment along the way will ever be wasted. If you buy the simplest system now, no future options are off the table.

Examples of the V-3 pump in use

The V-3 pump used in the R2-S is well proven. Nearly 200 of the Vapourtec V-3 pumps are in daily use in laboratories around the world. The V-3's capability to pump organometallic reagents continuously has been well documented in peer reviewed publications such as:

Continuous Flow-Processing of Organometallic Reagents Using an Advanced Peristaltic Pumping System and the Telescoped Flow Synthesis of (E/Z)-Tamoxifen:



<http://www.vapourtec.co.uk/products/eseriessystem/pumping/organometallic>

For more details visit www.vapourtec.com Tel: +44(0)1284 728659

Email: info@vapourtec.com