

# **New Applications of Flow Chemistry in Industry and Academia**

Dr Ryan Skilton MRSC  
Research Scientist, Vapourtec Ltd

[\*\*Ryan.Skilton@Vapourtec.com\*\*](mailto:Ryan.Skilton@Vapourtec.com)

## Vapourtec Key Facts

- Established 2003
- First Flow Chemistry System sold – May 2006
- Installed 370 systems worldwide
- Vapourtec only produce flow chemistry systems
- 250 peer-review publications using our systems to date




## Vapourtec Systems

- Vapourtec manufacture laboratory scale flow chemistry systems
  - R-Series (launched 2006); a highly specified modular system capable of either standalone operation or integration with other equipment to provide a versatile automated flow chemistry
  - E-Series (launched 2012); an easy to use entry level system platform ideal for new users and University teaching



# Where does Vapourtec fit in ?

***Flow principles can be applied to >14 orders of magnitude in scale, from earliest chemistry to late manufacturing***



50 mg (convenience)  
to 4000 g / day ("real chemistry")

## The Vapourtec R-Series (launched 2006)

- Highly specified modular system
- Capable of integrating with other equipment
- Pump options up to 200 bar
- Same reactor range as E-Series
- Either 2, 3 or 4 independent pumps
- 4 independently set reactors
- 270 systems installed worldwide
- Features in >200 peer reviewed publications





## The Vapourtec E-Series (launched Dec 2012)

- Entry-level system
- Ideal for academic groups
- Stand alone operation only
- V-3 pumps up to 10 bar pressure
- Same reactor range as R-Series
- Either 2 or 3 independent pumps
- 2 independent reactor positions
- Does not interface with other equipment
- 110 systems installed to date




## Vapourtec E-Series – The V-3 pump

- Pressure
  - 10 bar
- Compatibility
  - Strong acids and bases
  - Organometallics
  - Suspensions/slurries
  - Gases
- Smooth



## Vapourtec V-3 Pump – Advanced control

After




Pressure and flow rate of peristaltic V-3 pump showing the benefit of Vapourtec Adaptive Pump Rotor Control




## **Applications in Continuous Flow**

## Biphasic Reactions in Continuous Flow



**Laminar flow** or streamline flow in pipes (or tubes) occurs when a fluid flows in parallel layers, with no disruption between the layers. At low velocities, the fluid flows without lateral mixing, adjacent layers slide past one another like playing cards.


## Biphasic Reactions in Continuous Flow




**Static mixing inserts**

[Static Mixer Reactor](#)

## Biphasic Steven's Oxidation



One example of a reaction with a challenging mixing behaviour is an ester synthesis via a Steven's oxidation, which uses aqueous sodium hypochlorite to oxidise an aldehyde to an ester in the presence of an alcohol, allowing the direct synthesis of the ester from the aldehyde.






## Biphasic Steven's Oxidation


- Vapourtec Reactors for Rapid Mixing and **High Flow rate** pump module
- Greater than **99.9%** conversion at a **58 s** residence time in a **60 ml** reactor
- Over **4.5 kg/day** production with just the footprint of a Vapourtec R-Series

Reactor	Residence Time (s)	Conversion (%)	Throughput (gh <sup>-1</sup> )
10 ml Standard	58	64.6	22.4
20 ml Mixed	58	94.5	67.8
60 ml Mixed	58	99.9	203.3

## Biphasic Amide Formation



A method developed to permit amide formation without the need for 2 equivalents of amine – biphasic separation prevents reaction of base with the acid chloride



## Biphasic Amide Formation

- Greater than **95%** yield from a biphasic reaction
- **14.5 g/h** at 50 °C
- Useful synthetic strategy accelerated using intense mixing

Vapourtec wish to thank Dr Claudio Battilocchio and Prof. Steven Ley at the University of Cambridge for conducting the experimental work and providing additional information for this application note.

## **Solids in Flow?!**

- Common concern about blocking the flow path
- Difficult to achieve pressure – blocking pressure regulators




	TEMPERATURE (°C)															
	100	110	120	130	140	150	160	170	180	190	200	210	220	230	240	250
1 Butanol	0.5	0.8	1.1	1.5	2.1	2.8	3.7	4.7	6.1	7.7	9.5	11.7	14.2	17.1	20.4	24.1
1 Propanol	1.2	1.6	2.3	3.2	4.1	5.4	7.0	9.0	11.3	14.0	17.2	20.9	25.2	30.0	35.5	41.7
Acetic Acid	0.6	0.8	1.1	1.4	1.9	2.5	3.2	4.1	5.1	6.3	7.8	9.5	11.4	13.7	16.3	19.2
Acetone	3.7	4.8	6.1	7.6	9.5	11.6	14.1	16.9	20.1	23.8	27.9	32.5	37.6	43.2	49.4	56.2
Acetonitrile	2.0	2.7	3.7	4.8	6.3	8.1	10.3	13.0	16.3	20.2	24.8	30.2	36.6	44.0	52.5	62.3
Benzene	1.8	2.5	3.1	3.9	4.9	5.9	7.2	8.6	10.3	12.2	14.4	16.9	19.6	22.7	26.1	29.9
Carbon tetrachloride	1.9	2.5	3.2	4.0	5.0	6.1	7.4	8.9	10.5	12.4	14.6	16.9	19.6	22.4	25.6	29.0
Chloroform	3.1	4.0	5.0	6.3	7.8	9.6	11.7	14.1	16.8	19.9	23.4	27.3	31.7	36.6	42.0	47.9
Cyclohexane	1.7	2.2	2.9	3.6	4.5	5.5	6.7	8.1	9.7	11.5	13.5	15.7	18.3	21.0	24.1	27.5
DCM	5.9	7.5	9.4	11.7	14.3	17.4	20.9	24.9	29.4	34.5	40.2	46.5	53.5	61.2	69.6	78.7
Di ethyl ether	6.1	7.6	9.4	11.5	14.0	16.8	20.0	23.6	27.7	32.2	37.3	42.9	49.0	55.7	63.1	71.0
Diglyme	0.1	0.2	0.3	0.4	0.6	0.8	1.0	1.3	1.7	2.2	2.8	3.6	4.4	5.5	6.7	8.1
Dioxane	1.0	1.3	1.7	2.3	2.9	3.7	4.7	5.9	7.2	8.8	10.6	12.7	15.2	17.9	21.0	24.4
DME	1.2	1.6	2.0	2.5	3.0	3.7	4.5	5.4	6.4	7.5	8.8	10.2	11.8	13.5	15.4	17.5
DMF	0.2	0.3	0.4	0.6	0.8	1.0	1.3	1.7	2.1	2.6	3.2	3.9	4.7	5.6	6.6	7.8
DMSO	0.1	0.1	0.1	0.2	0.2	0.3	0.4	0.7	0.9	1.2	1.5	2.0	2.5	3.2	4.0	5.0
Ethanol	2.1	2.9	4.0	5.3	7.0	9.1	11.6	14.7	18.3	22.6	27.6	33.4	40.1	47.7	56.3	65.9
Ether	6.1	7.6	9.4	11.5	14.0	16.8	20.0	23.6	27.7	32.2	37.3	42.9	49.0	55.7	63.1	71.0
Ethyl Acetate	2.0	2.7	3.5	4.4	5.6	6.9	8.5	10.4	12.5	14.9	17.6	20.7	24.1	27.9	32.1	36.8
Formic Acid	1.0	1.3	1.7	2.2	2.9	3.6	4.5	5.6	6.9	8.3	10.0	11.9	14.1	16.6	19.4	22.5
Heptane	1.2	1.7	2.2	2.9	3.8	4.8	6.1	7.6	9.5	11.6	14.1	16.9	20.2	24.0	28.2	33.0
Hexane	2.5	3.1	4.0	5.0	6.1	7.4	9.0	10.7	12.7	14.9	17.4	20.2	23.2	26.5	30.1	34.1
IPA	2.0	2.8	3.8	5.1	6.7	8.7	11.0	13.9	17.2	21.2	25.7	30.9	36.8	43.5	51.0	59.4
MEK	1.9	2.4	3.1	4.0	4.9	6.1	7.5	9.0	10.8	12.8	15.1	17.6	20.4	23.4	26.8	30.5
MeOH	3.3	4.5	6.0	7.9	10.2	13.0	16.5	20.6	25.4	31.1	37.7	45.3	54.0	63.9	75.1	87.7
NMP	0.0	0.0	0.1	0.1	0.1	0.2	0.3	0.4	0.5	0.7	1.0	1.2	1.6	2.0	2.5	3.1
Pentane	5.9	7.3	9.0	10.9	13.0	15.5	18.3	21.4	24.8	28.5	32.7	37.2	42.0	47.3	52.9	59.0
p-Xylene	0.3	0.4	0.6	0.8	1.1	1.4	1.7	2.2	2.7	3.4	4.1	5.0	6.0	7.1	8.4	9.9
t Butyl Alcohol	1.9	2.7	3.6	4.8	6.3	8.1	10.2	12.7	15.6	19.0	22.9	27.3	32.3	37.8	44.0	50.8
THF	2.7	3.5	4.4	5.5	6.9	8.4	10.1	12.2	14.4	17.0	19.9	23.1	26.6	30.5	34.7	39.4
Toluene	0.7	1.0	1.3	1.7	2.2	2.7	3.4	4.2	5.1	6.2	7.5	8.9	10.5	12.3	14.4	16.7
Water	1.0	1.2	1.8	2.6	3.7	5.0	6.6	8.5	10.8	13.5	16.5	20.0	23.8	28.1	32.8	37.9


## Solids in Flow?!


- Common concern about blocking the flow path
- Difficult to achieve pressure – blocking pressure regulators
- Vapourtec have been using solids in slurries at pressure using V-3 pumps

## Palladium on Charcoal Slurries



## Palladium on Charcoal Slurries








## Palladium on Charcoal Slurries

- Over **6 g/h** of an API intermediate
- Pumping slurries of 5 and 10% palladium on charcoal **under pressure** continuously
- Heterogeneously catalysed transfer hydrogenation without catalyst scale limitations enabling straightforward scale-up
- **81%** isolated yield
- Use of the V-3 pump to control back pressure
- Versatile ability to optimise catalyst conditions


## Sodium Hydride as a Slurry



## Sodium Hydride as a Slurry

- Sodium hydride is one of the most versatile heterogeneous strong bases available in chemistry
- Many commercial forms of the reagent are cheaply available
- Commonly available 60% dispersion in mineral oil, a relatively stable and easily handled form
- Sodium hydride derived anions frequently lead to enhancements in yield and selectivity in a large variety of reactions


## Sodium Hydride as a Slurry



## Sodium Hydride as a Slurry

- Utilising sodium hydride as a slurry in THF for a laboratory scale continuous flow reaction
- Pumping the sodium hydride as a slurry and controlling gas formation by operation at a pressure of **7 bar**
- Truly continuous operation demonstrated during a **2-hour** synthesis of an alkyl carbazole heterocycle intermediate
- Fast reaction time of **4 minutes**
- 82 % isolated yield, 10.6 gh<sup>-1</sup> with greater than **98%** purity

## Singlet Oxygen Oxidation




- Oxygen is an ideal oxidising agent
- Clear potential dangers in its use at large scale
- Photochemical reactions can present scale-up challenges due to falling light penetration



## Features of the UV-150

- Fits both E and R-Series systems
- Multiple gram / hour scale-up
- High intensity UV light source
- User selectable UV power
- Light source wavelength filtering
- Temperature control - 40° C to 80° C
- Easily changed reactors
- Space saving compact design
- Interlocks ensure safe operation
- Optional spectrometer for real time monitoring of transmission spectra
- LED light sources are also available from 365 nm to 530 nm





## Singlet Oxygen Oxidation

- Pumping of reactive gases without the need for mass flow controllers
- Variable, automated back-pressure regulation
- Ease of scale-up with the UV-150
- Variable wavelength options – target specific chromophores or activation pathways
- **15.3 g/h**

## Industry and Academia

Vapourtec users include industrial ...

- Abbot, AstraZeneca, BP, GSK, J&J, Merck, Novartis, Pfizer, Roche, Sanofi-Aventis, UCB, Sigma Aldrich, Galderma, Novasep, Amgen, Takeda + 20 other Pharma, CROs and CMOs.

... and Academic

- Cambridge, MIT, Montreal, London (UCL), Nottingham, Cardiff Southampton, Toulouse, Bath, Paris SUD, FUNDP, Sheffield, Rennes, Milan, Kansas, CSIRO, Melbourne, Max Planck, + 10 other academic institutes.

The Vapourtec has been the system of choice for 250 peer reviewed publications to date.


370 Vapourtec systems installed worldwide.


<http://vapourtec.co.uk/publications>

## Pharmaceutical Industry

Janssen

Continuous Flow  $\alpha$ -Arylation of N,N-Dialkylhydrazones under Visible Light Photoredox Catalysis – Andrés A. Trabanco, 2017







## Pharmaceutical Industry

- Possible to obtain conversions of 89% in batch, with high selectivity
- On scale-up, conversion fell to **54%, and 50%** selectivity
- Translate batch optimal conditions into flow
- Scale-up of optimal gave comparable results in flow, permitting **91% conversion, and 75% isolated yield**

*“To our delight, the results were comparable when the reaction was scaled up to 15 mmol, and the product was isolated in 75% yield”*

## Academia

### Continuous Synthesis of Organozinc Halides Coupled to Negishi Reactions.


Organoboranes (Suzuki) are the most common choice in C-C bond formation and a wide range are **commercially available**

Organozinc (Negishi) are **commercially limited in range**


- superior reaction rates
- can be used when organoboranes are too unreactive
- can facilitate the coupling of  $sp^3$  carbons

However organozinc compounds can have issues with reproducibility and handling.

**If we could generate these intermediates and react them in-situ then we can increase the scope of C-C cross coupling diversity.**



## Academia



Long term stability assessment

- Column filled 12 g Zinc
- Initial trial used 5 g Zinc – 150 ml (0.5M) solution of benzyl zinc bromide (94%)
- Calculated turnover – 175
- Continuous output – 3.3 mmol h<sup>-1</sup>

**Approach is stable and robust enough to support larger scale chemistry**

# Academia


## Ethyl Lithiodiazoacetate: Extremely Unstable Intermediate Handled Efficiently in Flow



Diazo compounds are highly versatile, and ethyl diazoacetate (EDA) is a common diazo reagent, that can be used as a Nucleophile under base catalysis.

Aldehydes are sufficiently electrophilic to react with EDA directly, but other carbonyl electrophiles require the use of Organometallics, such as lithium diisopropylamide (LDA), generating ethyl lithiodiazoacetate.

## Academia





# Reactors



Standard PFA coil  
Ambient to 150 °C



High Temp coil  
Ambient to 250 °C



Cooled PFA coil  
Ambient to -70 °C



Standard Column  
Ambient to 150 °C



Cooled column  
Ambient to -40 °C



Gas/Liquid reactor  
Ambient to 150 °C



Micromixer reactor -40 °C to 150 °C



Photochemical reactor – 40 °C to 80 °C



Heated mixer reactor ambient to 150 °C

## Vapourtec Online

- Application notes
- Newsletter
- Recent Literature Publications
- Products

<https://www.vapourtec.com>

# **New Applications of Flow Chemistry in Industry and Academia**

Dr Ryan Skilton MRSC  
Research Scientist, Vapourtec Ltd

[\*\*Ryan.Skilton@Vapourtec.com\*\*](mailto:Ryan.Skilton@Vapourtec.com)