Catalytic processes under continuous flow-Current trends and applications

Manuel Nuño Chief Scientific Officer Vapourtec Ltd

manuel.nuno@vapourtec.com



Vapourtec Key Facts

- Established 2003
- First Flow Chemistry System sold May 2006
- Installed Over 500 R-Series & 230 E-Series systems
- More than 300 UV-150 photochemical reactors in use
- Vapourtec only produce flow chemistry systems
- Strong patent portfolio
- New product pipeline in development





Important features of flow chemistry

- Excellent Heat transfer surface area to volume ratio
- High rates of mass transfer small sizes
- Precise control of reaction parameters
- Low back-mixing = separation of products from reactants
- Only small quantities of hazardous materials "in-process"
- Reactive intermediates don't need to be isolated
- Zero head space = effective & safe at high pressure
- Run for longer to scale-up reactions

Flow chemistry – Basic set-up



Flow chemistry – Superheating of solvents

| | TEMPERATURE (°C) | | | | | | | | | | | | | | | |
|----------------------|------------------|-----|-----|------|------|------|------|------|------|------|------|------|------|------|------|------|
| | 100 | 110 | 120 | 130 | 140 | 150 | 160 | 170 | 180 | 190 | 200 | 210 | 220 | 230 | 240 | 250 |
| 1 Butanol | 0.5 | 0.8 | 1.1 | 1.5 | 2.1 | 2.8 | 3.7 | 4.7 | 6.1 | 7.7 | 9.5 | 11.7 | 14.2 | 17.1 | 20.4 | 24.1 |
| 1 Propanol | 1.2 | 1.6 | 2.3 | 3.2 | 4.1 | 5.4 | 7.0 | 9.0 | 11.3 | 14.0 | 17.2 | 20.9 | 25.2 | 30.0 | 35.5 | 41.7 |
| Acetic Acid | 0.6 | 0.8 | 1.1 | 1.4 | 1.9 | 2.5 | 3.2 | 4.1 | 5.1 | 6.3 | 7.8 | 9.5 | 11.4 | 13.7 | 16.3 | 19.2 |
| Acetone | 3.7 | 4.8 | 6.1 | 7.6 | 9.5 | 11.6 | 14.1 | 16.9 | 20.1 | 23.8 | 27.9 | 32.5 | 37.6 | 43.2 | 49.4 | 56.2 |
| Acetonitrile | 2.0 | 2.7 | 3.7 | 4.8 | 6.3 | 8.1 | 10.3 | 13.0 | 16.3 | 20.2 | 24.8 | 30.2 | 36.6 | 44.0 | 52.5 | 62.3 |
| Benzene | 1.8 | 2.5 | 3.1 | 3.9 | 4.9 | 5.9 | 7.2 | 8.6 | 10.3 | 12.2 | 14.4 | 16.9 | 19.6 | 22.7 | 26.1 | 29.9 |
| Carbon tetrachloride | 1.9 | 2.5 | 3.2 | 4.0 | 5.0 | 6.1 | 7.4 | 8.9 | 10.5 | 12.4 | 14.6 | 16.9 | 19.6 | 22.4 | 25.6 | 29.0 |
| Chloroform | 3.1 | 4.0 | 5.0 | 6.3 | 7.8 | 9.6 | 11.7 | 14.1 | 16.8 | 19.9 | 23.4 | 27.3 | 31.7 | 36.6 | 42.0 | 47.9 |
| Cyclohexane | 1.7 | 2.2 | 2.9 | 3.6 | 4.5 | 5.5 | 6.7 | 8.1 | 9.7 | 11.5 | 13.5 | 15.7 | 18.3 | 21.0 | 24.1 | 27.5 |
| DCM | 5.9 | 7.5 | 9.4 | 11.7 | 14.3 | 17.4 | 20.9 | 24.9 | 29.4 | 34.5 | 40.2 | 46.5 | 53.5 | 61.2 | 69.6 | 78.7 |
| Di ethyl ether | 6.1 | 7.6 | 9.4 | 11.5 | 14.0 | 16.8 | 20.0 | 23.6 | 27.7 | 32.2 | 37.3 | 42.9 | 49.0 | 55.7 | 63.1 | 71.0 |
| Diglyme | 0.1 | 0.2 | 0.3 | 0.4 | 0.6 | 0.8 | 1.0 | 1.3 | 1.7 | 2.2 | 2.8 | 3.6 | 4.4 | 5.5 | 6.7 | 8.1 |
| Dioxane | 1.0 | 1.3 | 1.7 | 2.3 | 2.9 | 3.7 | 4.7 | 5.9 | 7.2 | 8.8 | 10.6 | 12.7 | 15.2 | 17.9 | 21.0 | 24.4 |
| DME | 1.2 | 1.6 | 2.0 | 2.5 | 3.0 | 3.7 | 4.5 | 5.4 | 6.4 | 7.5 | 8.8 | 10.2 | 11.8 | 13.5 | 15.4 | 17.5 |
| DMF | 0.2 | 0.3 | 0.4 | 0.6 | 0.8 | 1.0 | 1.3 | 1.7 | 2.1 | 2.6 | 3.2 | 3.9 | 4.7 | 5.6 | 6.6 | 7.8 |
| DMSO | 0.1 | 0.1 | 0.1 | 0.2 | 0.2 | 0.3 | 0.4 | 0.7 | 0.9 | 1.2 | 1.5 | 2.0 | 2.5 | 3.2 | 4.0 | 5.0 |
| Ethanol | 2.1 | 2.9 | 4.0 | 5.3 | 7.0 | 9.1 | 11.6 | 14.7 | 18.3 | 22.6 | 27.6 | 33.4 | 40.1 | 47.7 | 56.3 | 65.9 |
| Ether | 6.1 | 7.6 | 9.4 | 11.5 | 14.0 | 16.8 | 20.0 | 23.6 | 27.7 | 32.2 | 37.3 | 42.9 | 49.0 | 55.7 | 63.1 | 71.0 |
| Ethyl Acetate | 2.0 | 2.7 | 3.5 | 4.4 | 5.6 | 6.9 | 8.5 | 10.4 | 12.5 | 14.9 | 17.6 | 20.7 | 24.1 | 27.9 | 32.1 | 36.8 |
| Formic Acid | 1.0 | 1.3 | 1.7 | 2.2 | 2.9 | 3.6 | 4.5 | 5.6 | 6.9 | 8.3 | 10.0 | 11.9 | 14.1 | 16.6 | 19.4 | 22.5 |
| Heptane | 1.2 | 1.7 | 2.2 | 2.9 | 3.8 | 4.8 | 6.1 | 7.6 | 9.5 | 11.6 | 14.1 | 16.9 | 20.2 | 24.0 | 28.2 | 33.0 |
| Hexane | 2.5 | 3.1 | 4.0 | 5.0 | 6.1 | 7.4 | 9.0 | 10.7 | 12.7 | 14.9 | 17.4 | 20.2 | 23.2 | 26.5 | 30.1 | 34.1 |
| IPA | 2.0 | 2.8 | 3.8 | 5.1 | 6.7 | 8.7 | 11.0 | 13.9 | 17.2 | 21.2 | 25.7 | 30.9 | 36.8 | 43.5 | 51.0 | 59.4 |
| MEK | 1.9 | 2.4 | 3.1 | 4.0 | 4.9 | 6.1 | 7.5 | 9.0 | 10.8 | 12.8 | 15.1 | 17.6 | 20.4 | 23.4 | 26.8 | 30.5 |
| MeOH | 3.3 | 4.5 | 6.0 | 7.9 | 10.2 | 13.0 | 16.5 | 20.6 | 25.4 | 31.1 | 37.7 | 45.3 | 54.0 | 63.9 | 75.1 | 87.7 |
| NMP | 0.0 | 0.0 | 0.1 | 0.1 | 0.1 | 0.2 | 0.3 | 0.4 | 0.5 | 0.7 | 1.0 | 1.2 | 1.6 | 2.0 | 2.5 | 3.1 |
| Pentane | 5.9 | 7.3 | 9.0 | 10.9 | 13.0 | 15.5 | 18.3 | 21.4 | 24.8 | 28.5 | 32.7 | 37.2 | 42.0 | 47.3 | 52.9 | 59.0 |
| p-Xylene | 0.3 | 0.4 | 0.6 | 0.8 | 1.1 | 1.4 | 1.7 | 2.2 | 2.7 | 3.4 | 4.1 | 5.0 | 6.0 | 7.1 | 8.4 | 9.9 |
| t Butyl Alcohol | 1.9 | 2.7 | 3.6 | 4.8 | 6.3 | 8.1 | 10.2 | 12.7 | 15.6 | 19.0 | 22.9 | 27.3 | 32.3 | 37.8 | 44.0 | 50.8 |
| ТНЕ | 2.7 | 3.5 | 4.4 | 5.5 | 6.9 | 8.4 | 10.1 | 12.2 | 14.4 | 17.0 | 19.9 | 23.1 | 26.6 | 30.5 | 34.7 | 39.4 |
| Toluene | 0.7 | 1.0 | 1.3 | 1.7 | 2.2 | 2.7 | 3.4 | 4.2 | 5.1 | 6.2 | 7.5 | 8.9 | 10.5 | 12.3 | 14.4 | 16.7 |
| Water | 1.0 | 1.2 | 1.8 | 2.6 | 3.7 | 5.0 | 6.6 | 8.5 | 10.8 | 13.5 | 16.5 | 20.0 | 23.8 | 28.1 | 32.8 | 37.9 |

Flow chemistry – Multi-step reactions



- Intermediate may be
 - Unstable
 - Toxic
 - Air or moisture sensitive
- Used immediately
 - Minimal inventory

Where our Customers apply flow chemistry

- Highly exothermic reactions
- Reactions involving unstable (or toxic) intermediates
- Very rapid reactions (Rt < 1 min)
- Reactions requiring superheating
- One or more very volatile reagents or dissolved gases
- Multi-phase: liquid/liquid, liquid/solid or gas/liquid
- Reactions requiring better selectivity
- Continuous flow photochemistry
- Continuous flow electrochemistry
- Easy scaleup route once reaction is optimized

Catalysis in flow – How to design a reaction

To translate batch catalytic reactions in continuous flow, three questions must be answered:

| • | Physical form of catalyst? | It will define how the catalyst is handled |
|---|---|---|
| • | Physical properties of reagents? | It will define how to handle them and reaction parameters |
| • | Thermally or non-thermally mediated reaction? | It will define the type of reactor needed |

Catalysis in flow – Physical form of catalyst

 Homogeneous catalysts – Same phase as reagents, usually liquid

 Heterogeneous catalyst –Different phase as reagent, usually solid





Catalysis in flow – Homogeneous catalysis

- Catalyst and reagent premixed in the same solution
 - Simple approach
 - Fixed stoichiometric ratios per solution



ourtec

- Catalyst and reagent prepared in different solutions
 - Multi-channel approach
 - Tunable stoichiometric ratios with the same solutions



urtec

Catalysis in flow – Homogeneous catalysis

Decarboxylative cross coupling – Substrate scope*





| Batch | 7 % | 37 % |
|-------|------|------|
| Flow | 71 % | 0 % |



70% 60%



NO₂

71%









*Vapourtec Application Note 27 : Decarboxylative Cross-Couplings with a Soluble Catalyst System

Catalysis in flow – Heterogeneous catalysis

- Fixed bed reactor Solid catalyst packed in a column, reagents flow through (*i.e. immobilized enzymes, Raney-Ni catalyst, etc.*)
- Solid catalyst pumped as slurry Some catalysts can be pumped as slurry (*i.e. Pd/C, TiO₂, etc.*)
- Catalytic reactors The reactor itself is the catalyst (*i.e. Copper reactor*)







Catalysis in flow – Heterogeneous catalysis

Suzuki-Miyaura coupling*







*Vapourtec Application Note 49 – Suzuki Coupling with SiliaCat DPP-Pd Heterogeneous Catalyst

Catalysis in flow – Heterogeneous catalysis

Suzuki-Miyaura coupling*



*Vapourtec Application Note 49 – Suzuki Coupling with SiliaCat DPP-Pd Heterogeneous Catalyst

Catalysis in flow – Heterogeneous catalysis

Palladium on Charcoal Slurries with H₂ Gas*





- Pd/C 5 and 10 %
- Up 100 mg/ml
- Pressure with peristaltic BPR

- H₂ using same pump type
- 16 experiments in under 12 h
- Selectivity

Catalysis in flow – Heterogeneous catalysis

Palladium on Charcoal Slurries with H₂ Gas*

- Reagents prepared in a single reservoir
- Pd/C suspended in EtOAc on a stirrer plate prevents settling
- Peristaltic V-3 is able to pump the slurry directly
- T-piece arranged so that the slurry doesn't change direction
- Peristaltic BPR makes it possible to use the slurry at pressures up to 10 bar



vapourtec

*Vapourtec Application Note 54 – Selective hydrogenation of O-benzyl vanillin using hydrogen gas and a palladium on charcoal slurry

Catalysis in flow – Heterogeneous catalysis

urtec

Palladium on Charcoal Slurries with H₂ Gas*

- Pumping slurries of up 100 mg/ml 5 and 10% palladium on charcoal under pressure continuously
- 3-phases in flow
- Avoids pressure limitations of packed beds
- 81% isolated yield selective between products
- Use of the V-3 pump to control back pressure
- Versatile ability to optimise catalyst conditions
- Over **6** g/h of an API intermediate from a 1 mm bore reactor



*Vapourtec Application Note 54 – Selective hydrogenation of O-benzyl vanillin using hydrogen gas and a palladium on charcoal slurry

Catalysis in flow – Novel reaction pathways

What happens when the reaction is not thermally mediated?

- Photocatalysts need of photons to promote reactions
- Electrocatalysts need of electrons to promote reactions





ourtec

- These reactions are traditionally difficult in batch due to:
 - Non-homogeneous radiation field
 - Hot-spots

```
• ...
```

Catalysis in flow – Photocatalysis

Different photocatalysts require of different wavelengths

- TiO₂, TBADT ~365 nm
- 4CzIPN ~ 420 nm
- [Ir(dFCF3ppy)₂(5,5'-dCF₃bpy)](PF₆)~450 nm



Low pressure mercury - 3 options:

- 254 nm
- 310 nm
- 370 nm

Medium pressure mercury: 220 nm to 600 nm filters to select desired wavelength

LEDs:

a range of precise wavelengths 365-700 nm 60 or 150 W input power







Catalysis in flow – Photocatalysis

Photocatalytic synthesis of γ-lactams*

- Simple and direct synthesis of γ -lactams and α -tertiary amine derivative
- (HAT) catalyst using a cheap organic photocatalyst (4CzIPN) in combination with azide ion Other N-Functionalisations of Flow-Generated y-Amino Esters Fmoc CO₂Me NH CO₂Me CO₂Me NH Bu₄t (10 m 36, 63% 37, 58% 38, 52% (1.0 equiv, 0.3 M) auiv) [N-protection] [amidation] [amidation] CO₂Me 4Cz 2h 2a (1 m y-lactam 4a CO₂t-Bu (1.0 equiv, 0.3 M) NH CO₂Me F 39.50% 40, 41% [sulfonamidation] [reductive amination]

*1. Vapourtec Application Note 68: Photocatalytic Synthesis of γ-Lactams and α-Tertiary Amine Derivatives in Continuous Flow *2. A. S. H. Ryder et al. Angew. Chem. Int. Ed., 2020

Catalysis in flow – Photocatalysis

Photocatalytic synthesis of γ-lactams – Effect of photonic loading*



*1. Vapourtec Application Note 68: Photocatalytic Synthesis of γ-Lactams and α-Tertiary Amine Derivatives in Continuous Flow *2. A. S. H. Ryder et al. Angew. Chem. Int. Ed., 2020

Catalysis in flow – Electrocatalysis

- Integrated or stand-alone versions
- Working temperatures: from -10 °C to 100 °C
- Limit pressure up to 5 bar, allowing to work above solvent's boiling point and with gas mixtures
- Vapourtec supplies 20 different electrodes







Catalysis in flow – Electrocatalysis

Electrochemistry – Basic principles

- Reactions carried on an electrolytic cell
- Flow of electrons drives non-spontaneous redox reactions
- Working in batch held back electrosynthesis...



Non-homogeneous electric field

Need of electrolyte

Hot-spots on solution media

Low efficiency and/or difficult purification

Catalysis in flow – Electrocatalysis

New opportunities in flow

- By flowing reagents through a microreactor you can overcome batch issues
 - Homogeneous electric field
 - Electrolyte is now an option
 - Homogeneous temperature
- Why to switch from organic oxidants?
 - Non-hazardous reagents
 - Selectivity based on current



Catalysis in flow – Electrocatalysis

Electrochemical Csp2-Csp3 cross coupling of organic halides*





| | Reagents | lon electrochemica | I | Collection of products via liquid handler | Entry | Temperature | Flow rate | Rt | Ligand | SM | тм |
|--|----------------------|-----------------------|-------------|---|-------|-------------|-------------|--------|--------|------|------|
| | & Catalyst | reactor | Passive BPR | | 1 | 30 °C | 0.06 ml/min | 10 min | Dtbbpy | 38 % | 18 % |
| | | | | | 2 | 30 °C | 0.03 ml/min | 20 min | Dtbbpy | 5 % | 40 % |
| | Ŷ | | • | | 3 | 50 °C | 0.06 ml/min | 10 min | Dtbbpy | 9 % | 38 % |
| | | | | | 4 | 50 °C | 0.03 ml/min | 20 min | Dtbbpy | 0 % | 50 % |
| | $\square \checkmark$ | I I ISI I I | | | 5 | 50 °C | 0.06 ml/min | 10 min | L1 | 5 % | 55 % |
| | | | | | 6 | 50 °C | 0.03 ml/min | 20 min | L1 | 0 % | 81 % |

*Application Note 63: Electrochemical pathway for cross coupling of organic halides - Csp2-Csp3 bonding

Thank you for your time today

- Website
 <u>www.vapourtec.com</u>
- Contact manuel.nuno@vapourtec.com
- Publications
 <u>www.vapourtec.co.uk/publications</u>
- Newsletter
 - Latest http://www.vapourtec.co.uk/newsletter/latestissue
 - Signup http://www.vapourtec.co.uk/newsletter/signup

